

DigiPREP HT

SULFURIC ACID (H₂SO₄) DIGESTION & DISTILLATION SYSTEMS




SAMPLE DIGESTIONS

The *Digi*PREP HT series of block digestion systems offers the latest in graphite block technology in four different packages. Graphite blocks are coated with silicon carbide (SiC) which resists aggressive, corrosive attack for guaranteed long life in harsh laboratory environments. The systems have a temperature range from ambient to 450 °C and require a *Digi*PREP controller.

DIGESTION MADE EASY...

The *Digi*PREP HT series (system(s) require the use of a controller) is available in four models - two models employing 100 ml tubes and two models with 250 ml tubes. An all-Teflon® manifold, with its unique Teflon® tube adaptors, removes corrosive acid fumes while promoting gentle refluxing. With its two-tier, dual zone tube rack, benchtop space is reduced by 50% allowing more room for other lab instruments. *Digi*PREP HT systems are ideal for Kjeldahl/TKN, WEEE RoHS, and for any application using sulfuric acid (H₂SO₄).



BUILD A COMPLETE SYSTEM

Description	DigiPREP HT 100-10	DigiPREP HT 100-40	DigiPREP HT 250-10	DigiPREP HT 250-20
Catalog Number	010-520-435	010-520-402	010-520-445	010-520-422
Capacity	10-100 ml	40-100 ml	10-250 ml	20-250 ml
Power	3000 Watt	3362 Watt	3000 Watt	3362 Watt
Size (in.)	17.75W x 12.5D x 27H	17W x 17.5D x 27H	17.75W x 12.5D x 27H	17W x 17.5D x 27H
Weight (lb/kg)	23 / 10.5	49 / 22	23 / 10.5	46 / 21

TOTAL CONTROL AT YOUR FINGERTIPS

The Color Touch Screen Controller (010-500-250)

is recommended for all multi-step, multi-temperature digestions such as Kjeldahl/TKN and WEEE RoHS digestions.

Color Touch Screen Controller with USB Port

(010-500-275). The USB port is conveniently located at the back of the controller and enables users to save detailed reports for laboratory audits. The reports are saved onto a USB stick and include the following data: date and time of the sample run, instrument serial number, temperature log in 15 second intervals and safety messages. The controller is perfect for research and method development.



Features

Allows for block temperature feed back or sample temperature feed back using an optional *DigiPROBE*

Automatic shut off at end of run

Real time graphical representation of heating progress

Displays temperature to 0.1 °C

Displays time to 0.1 seconds

10 languages included: English, French, German, Dutch, Italian, Portuguese, Spanish, Japanese, Chinese and Korean

Features

Temperature Set-Point Programming to 0.1 °C

Block sensor or *DigiPROBE* can be calibrated with up to 5 points



KeyPad Controller (010-500-220) is recommended for simpler methods where a single set-point temperature is required. Display shows current block temperature along with current time remaining for the digestion. Built-in audible and visual alarms signal the end of a segment. User friendly, one-touch buttons allow for fast resetting of temperature or alarm times.

The PC Classic Controller (010-500-280) and Go2Temp® software allow a user to control one or multiple *Digi*PREP blocks remotely. Access to advanced software features of the Go2Temp® is provided. Comes complete with communication capabilities (RS485 to USB).

System	Catalog Number	Qty	Component Description
<i>Digi</i> PREP HT 100-40 (230V) 010-520-402	010-520-010	1	<i>Digi</i> PREP HT 100-40 Insert Rack
	010-520-100	1	Exhaust Manifold Rack, 100 ml
Models for 100 ml digestion tubes	010-520-205	1	DigiPREP HT 100 (230V) block with Hi-rise
	010-520-300	1	Glass Water Jet
	010-520-304	300 gr	Boiling Stones
Disciplination of the second	010-520-230	1	DigiPREP HT 250-20 Insert Rack
DigiPREP HT 250-20 (230V) 010-520-422	010-520-225	1	DigiPREP HT 250 (230V) block with Hi-rise
Models for 250 ml digestion	010-520-250	1	Exhaust Manifold Rack, 250 ml
tubes	010-520-300	1	Glass Water Jet
	010-520-304	300 gr	Boiling Stones

System Catalog Number		Qty	Component Description	
	010-520-247	1	DigiPREP HT 100-10 Insert Rack	
DigiPREP HT	010-520-253	1	Teflon® Manifold HT 100-10	
100-10 (230V) 010-520-435	010-520-300	1	Glass Water Jet	
010 020 400	010-520-304	300 gr	Boiling Stones	
	010-520-211	1	DigiPREP HT 250-10 Insert Rack	
DigiPREP HT	010-520-215	1	Teflon® Manifold HT250-10	
250-10 (230V) 010-520-445	010-520-300	1	Glass Water Jet	
	010-520-304	300 gr	Boiling Stones	

ACCESSORIES

Description	Qty	Catalog No.
HT 100 Digestion Tubes, uncalibrated	1	010-520-011
HT 100 Digestion Tubes, calibrated * at the 25 ml and 50 ml mark	1	010-520-111
HT 100 Digestion Tubes, Quartz, uncalibrated	1	010-520-085
Glass Water Jet, Extra	1	010-520-300
Boiling Stones, Extra	300 gr	010-520-304
Ball Glass Condenser 100 ml	1	010-520-318

Description	Qty	Catalog No.
HT 250 Digestion Tubes, uncalibrated	1	010-520-031
HT 250 Digestion Tubes, calibrated * at the 50 ml and 100 ml mark	1	010-520-131
HT 250 Digestion Tubes, Quartz, uncalibrated	1	010-520-087
Ball Glass Condenser 250 ml	1	010-520-312
Heat deflector Shields for all racks	2	010-520-032

^{*} Custom tube calibrations available upon request.

UNIVERSAL FUME SCRUBBER (010-520-160)

The Universal Fume Scrubber is used to trap and neutralize fumes from the DigiPREP HT or from a DigiPREP MS with a DigiVAC Hood. Acid fumes travel through a series of bottles where the fumes are condensed, neutralized and scrubbed. An adapter is included to connect the system to either the HT Digestion manifold or the *Digi*VAC system.



DigiPREP distillation systems range from manual to fully automatic. Systems are programmed and operated using a single, turn-and-press knob. All systems can be used for the steam distillation of Kjeldahl ammonia solutions. In the Kjeldahl application, the additions of NaOH, H₂O, and H₂BO₂ may be automated. Samples can be distilled in 3 to 5 mintues. Variable steam power from 40% to 99% is offered for alcohol and phenol analysis in most models. Automatic drainage of hot sample waste provides safety and convenience. An instrument is available to suit every budget. Please review the chart below to determine model of choice.



pH TITRATOR (010-035-050)

Accurate dosing technology and high-precision glass cylinders ensure correct titration results. An easy to read display shows the quantity of the titrant used; delivers with 0.01 ml accuracy. May be connected to *DigiPREP* 500 distillation system to automate the final stage of Kjeldahl analysis.



	Model 100	Model 200	Model 300	Model 400	Model 500
Catalog Number	010-035-001	010-035-002	010-035-003	010-035-004	010-035-005
Automatic Addition of H ₂ O	No	Yes	Yes	Yes	Yes
Automatic Addition of NaOH	Yes	Yes	Yes	Yes	Yes
Automatic Addition of H ₃ BO ₃	No	No	No	Yes	Yes
Programmable Reaction Time	Yes	Yes	Yes	Yes	Yes
Programmable Distillation Time	Yes	Yes	Yes	Yes	Yes
Automatic Steam Generation	Yes	Yes	Yes	Yes	Yes
Variable Steam Generating capacity (40%-100%)	Yes	Yes	Yes	Yes	Yes
Automatic Removal of Sample Waste by Suction	No	No	Yes	Yes	Yes
Number of Programs	1	1	10	99	99

SERVICE DEPARTMENT

Our friendly Service Department will gladly assist you with technical questions pertaining to *DigiPREP* products. With over 30 years in servicing instruments, our staff will review and resolve issues quickly.

www.scpscience.com | sales@scpscience.com











CANADA / USA

Tel.: +1 (800) 361-6820 Fax: +1 (800) 253-5549 FRANCE

Tel.: +33 (0) 1 69 18 71 17 Fax: +33 (0) 1 60 92 05 67 **GERMANY**

Tel.: +49 (0) 8342-89560-61 Fax: +49 (0) 8342-89560-69 **INTERNATIONAL**

Tel.: +1 (514) 457-0701 Fax: +1 (514) 457-4499